

“Bottle Adapter“

Declaration according to article 22, 2017/745/EC

Technologie Institut Medizin GmbH (TIM), August - Thyssen - Str. 30, D-56070 Koblenz, is the legal manufacturer of the product “bottle adapter“ with models BA-ISO-K, BA-SEV-K, BA-SEV-Q, BA-SEV-SV, BA-SEV-PCC, BA-DES-S.

The "bottle adapter" (BA) product is used in the intended use to connect pharmaceutical containers with liquid volatile anesthetic agent and an anesthetic vaporizer to fill the anesthetic vaporizer with the liquid anesthetic agent contained in the bottle. The bottle adapters are coded specifically for the anesthetic gas and may only be used with the appropriately coded pharmaceutical containers.

List of anesthetic vaporizers compatible with “bottle adapter“:

Manufacturer Anesthetic vaporizer	Model Anesthetic vaporizer	Modell Bottle Adapter	Restrictions to be observed
Technologie Institut Medizin GmbH (TIM)	MIRUS Controller ISO	BA-ISO-K	None
Technologie Institut Medizin GmbH (TIM)	MIRUS Controller SEVO	BA-SEV-K BA-SEV-Q BA-SEV-SV BA-SEV-PCC	None
Technologie Institut Medizin GmbH (TIM)	MIRUS Controller DES	BA-DES-S	None

We, Technologie Institut Medizin GmbH (TIM), declare, that the listed combination of anesthetic vaporizer and “bottle adapter” BA was tested for mutual compatibility of the products for safety and their intended use according to the manufacturer's specifications and that these listed combinations can be used when any given restrictions and safety precautions in the instruction for use of the “bottle adapter” and the anesthetic vaporizer are observed.

Koblenz, den 21.05.2021



Prof. Ing. Thomas P.Kriesmer

://end of document

