

Service Bulletin

IMPORTANT INFORMATION – PLEASE READ AND KEEP



Change in Safety Check - STK

Document:	021_Change in STK_SB_03_A-00_EN
Date:	2020-05-25
Author:	Stephan Lang
Related version:	<input checked="" type="checkbox"/> All
Distribute to:	<input type="checkbox"/> Internal Use Only (Confidential) <input checked="" type="checkbox"/> Medcaptain <input checked="" type="checkbox"/> Dealer <input checked="" type="checkbox"/> End-user

STK adjustment from June 1st 2020

In order to cope with the increased number of necessary STKs without major time gaps and thus to maintain the operation of the MIRUS in operation continuously, an STK adjustment is required which will be **valid from June 1st 2020**. The valid protocol is version **C.00**. This can be found at <https://tim-gmbh.de/en/mirus-docs/>.

STK 12M: All necessary maintenance measures (A.-D.) of the twelve-month safety check can be carried out on site by any certified medical or service technician using common standard equipment and provided with his own seal/ stamp.

STK 24M: Every 24 months the major maintenance (A.-G.) must take place. This requires special test equipment and spare parts which can be purchased from TIM.



Of course, TIM is happy to continue to take on both the legally required STK 12M and the STK 24M. Here the known conditions apply.

//:End of the document

021_Change in STK_SB_03_A-00_EN

