

MIRUS™ & Bottle Adapter – Pharmaceutical list

Compatibility list of pharmaceutical manufacturers

Technologie Institut Medizin GmbH (TIM), August - Thyssen - Str. 30, D-56070 Koblenz, is the legal manufacturer of the products “MIRUS™” and “Bottle Adapter” (BA). The product “Bottle Adapter” (BA) is intended to connect pharmaceutical containers, containing liquid volatile anesthetic agent, and an anesthetic vaporizer (MIRUS) to fill the MIRUS™ with the liquid anesthetic agent contained in the bottle. The BA are key coded specifically for one anesthetic agent and may only be used with the appropriate key coded pharmaceutical containers.

Compatibility list of pharmaceutical manufacturers and substances to fill „MIRUS™“:

Pharmaceutical Manufacturer	Pharmaceutical substance	Modell Bottle Adapter	Market name
AbbVie	Isoflurane	BA-ISO-K	Forane
Baxter	Isoflurane	BA-ISO-K	Forane, Aerrane
Piramal Critical Care	Isoflurane	BA-ISO-K	Isoflurane, USP
AbbVie	Sevoflurane	BA-SEV-K	Ultane
Baxter	Sevoflurane	BA-SEV-K	Sevoflurane, USP
Baxter	Sevoflurane	BA-SEV-SV	Sevoflurane, USP (with integrated Safe-T-Seal valve)
Piramal Critical Care	Sevoflurane	BA-SEV-K	Sevoflurane, USP
AbbVie	Sevoflurane	BA-SEV-Q	Sevorane
Piramal Critical Care	Sevoflurane	BA-SEV-PCC	Sevoflurane Piramal 100% (with integrated safety valve)
Baxter	Desflurane	BA-DES-S	Suprane

We, Technologie Institut Medizin GmbH (TIM), declare, that the listed combination of the pharmaceutical manufacturer with the named substances and “Bottle Adapter” BA and “MIRUS™” was tested for compatibility and safe operation. All listed volatile substances can be used without hesitation in combination with the corresponding BA for safe filling of the MIRUS™ anesthesia system.

Koblenz, 2022-02-21



Prof. Ing. Thomas P.Kriesmer



//: End of document.