







ORS Clean Air





Passive patient gas disposal



Safety on a clear path




Compatibility Matrix




Model Clean Air	Ventilator		Adapter Clean Air		Image (e.g.)
	Manufacturer	Model	Art.-No.	Specifications	
ORS-CA-01 	Dräger Medical GmbH	Evita 2 Evita 2 dura Evita 2 cap Evita 4 Evita XL Savina 300	ORS-VC-04	22 mm F (to ventilator) / 30 mm M (to Clean air) Length: 75 cm	 Savina 300
	Dräger Medical GmbH	Evita Infinity V 300 Evita Infinity V 500 Evita Infinity V 600 Evita Infinity V 800	ORS-VC-09	22 mm F (to ventilator) / 30 mm M (to Clean air) Length: 37 cm	 Evita V 500
	GE Healthcare Inc.	Engström Carestation Carescape R860	ORS-VC-02	30 mm F (to ventilator) / 30 mm M (to Clean Air) Length: 75 cm	
	Hamilton Medical AG	C2 C3 C3S C6 C6S	ORS-VC-11	26 mm F (to ventilator) / 30 mm M (to Clean Air) Length: 75 cm	NOTE:  ORS-VC-11 is not to connect airtight to the cone of the exp. valve exit of the Hamilton ventilator but fits into the gap between the exp. valve exit

Model Clean Air	Ventilator		Adapter Clean Air		Image (e.g.)
	Manufacturer	Model	Art.-No.	Specifications	
ORS-CA-01 	Hamilton Medical AG	C1	ORS-VC-12	23 mm M (to ventilator) / 30 mm M (to Clean Air) Length: 75 cm	
	Maquet Critical Care AB	Servo-s	ORS-VC-08	30 mm F (to ventilator) / 30 mm M (to Clean air) Length: 37 cm	
	Medtronic (Nelcor, Puritan Bennett) Inc.	PB 840	ORS-VC-06	26 mm M (to ventilator) / 30 mm M (to Clean Air) Length: 75 cm	
	Fritz Stephan GmbH	EVEin	ORS-VC-18	22 mm M (to ventilator) / 30 mm M (to Clean Air) Length: 150 cm	

Model Clean Air	Ventilator		Adapter Clean Air		Image (e.g.)
	Manufacturer	Model	Art.-No.	Specifications	
<p>ORS-CA-02</p> 	Maquet Critical Care AB	Servo-i Servo-u	ORS-VC-08	30 mm F (to ventilator) / 30 mm M (to Clean air) Length: 37 cm	 <p>Servo-i</p>

Model Clean Air	Ventilator		Adapter Clean Air		Image (e.g.)
	Manufacturer	Model	Art.-No.	Specifications	
ORS-CA-04 	Dräger Medical GmbH	Evita 4 Evita XL Savina 300	ORS-VC-09	22 mm F (to ventilator) / 30 mm M (to Clean air) Length: 37 cm	 Evita XL
	Hamilton Medical AG	G5 and S1 with Expiration valve (reusable) G5 and S1 with Expiration valve (reusable) and with humidifier H 900	ORS-VC-10	30 mm F (to ventilator) / 30 mm M (to Clean air) Length: 37 cm	
	Hamilton Medical AG	C6 C6S	ORS-VC-16	26 mm F (to ventilator) / 30 mm M (to Clean air) Length: 37 cm	
	Löwenstein Medical Innovation GmbH & Co. KG	Elisa 600 Elisa 800	ORS-VC-17	27 mm F (to ventilator) / 30 mm M (to Clean air) Length: 37 cm	

Model Clean Air	AWS		Adapter Clean Air	General Note
	Manufacturer	Model	Art.-No.	
ORS-CA-01 	Dräger Medical GmbH	Julian Primus Cicero Cato Fabius Fabius Trio	ORS-VC-13 [Length: 75 cm]	The ORS-01 can be adapted to further AWS devices if the AWS has an AGSS port
	GE Healthcare Inc.	Aestiva Avance Aisys	ORS-VC-14 [Length: 113 cm] ORS-VC-15 [Length: 150 cm]	



Technologie Institut Medizin GmbH (TIM)
August-Thyssen-Str. 30
56070 Koblenz
Germany
Tel.: +49 261 899 689 00
Fax.: +49 261 899 689 09
www.tim-gmbh.de
info@tim-gmbh.de



Specifications are subjects to change without prior notice.
ORS Clean Air Compatibility Matrix ©2022 TIM GmbH.
All rights reserved. Rev. E-00 of 06/22.
TIM are a trademark of TIM GmbH.



Visit us
www.ors-clean-air.de

